

16 JUN 21 BRN  
Batch 1

Carey Reservoir 2021 and a few 2020 returns  
Solubles

p1 of 2

# RACK MAP--LACHAT FIA

1	spiked blank	16	F 24 MAY 21 6.2	31	F 14 MAY 21 3.8	46	B 08 MAR 21 6.0	61	F 08 FEB 21 8.0	76	F 26 FEB 21 5.0
2	"	17	F 24 MAY 21 wet	32	B 05 APR 21 11.0	47	F 16 APR 21 9.0	62	C 19 APR 21 11.0	77	B 05 APR 21 3.0
3	"	18	F 08 FEB 21 3.8	33	B 03 MAY 21 0.1	48	F 24 MAY 21 1.6	63	F 30 SEP 20 9.0	78	F 08 FEB 21 9.0
4	F 14 MAY 21 1.6	19	F 24 MAY 21 9.0	34	F 16 APR 21 8.0	49	F 22 MAR 21 0.1	64	F 05 OCT 20 wet	79	F 26 FEB 21 0.1
5	B 03 MAY 21 9.0	20	B 03 MAY 21 11.0	35	F 06 MAY 21 5.0	50	C 11 MAY 21 ceto.	65	B 20 AUG 20 10.0	80	B 14 MAY 21 11.0
6	F 16 APR 21 1.6	21	B 03 MAY 21 6.0	36	F 08 FEB 21 5.0	51	F 24 MAY 21 8.0	66	C 22 JUL 20 21.0	81	F 06 MAY 21 6.2
7	F 22 MAR 21 8.0	22	F 16 APR 21 6.2	37	F 24 MAY 21 0.1	52	F 22 MAR 21 3.8	67	F 29 JUN 20 0.1	82	B 14 MAY 21 3.0
8	B 05 APR 21 6.0	23	F 16 APR 21 weir	38	F 14 MAY 21 0.1	53	F 16 APR 21 5.0	68	F 02 SEP 20 1.6	83	F 24 MAY 21 3.8
9	F 22 MAR 21 9.0	24	C 19 APR 21 20.0	39	F 26 FEB 21 8.0	54	F 24 MAY 21 5.0	69	F 02 DEC 20 wet	84	F 08 FEB 21 1.6
10	F 26 FEB 21 1.6	25	C 19 APR 21 mfc 300	40	F 14 MAY 21 5.0	55	F 26 FEB 21 6.2	70	F 09 NOV 20 9.0	85	F 26 APR 21 wet 200
11		26		41		56		71		86	
12		27		42		57		72		87	
13		28		43		58		73		88	
14	Duplicate	29	Duplicate	44	Duplicate	59	Duplicate	74	Duplicate	89	Duplicate
15	Spike	30	Spike	45	Spike	60	Spike	75	Spike	90	Spike

Remember to do a duplicate, a check standard, and a blank in each set of 10 samples.

Spiking Info:

4mls sample + 80 µl 1000 ppb

Sample Info:

last std in middle batch  
had no tube - drew air

16 JUN 21

Carey Reservoir  
Solubles

cont

P2 of 2

## RACK MAP--LACHAT FIA

1	F 26 FEB 21 wet	16	F 16 APR 21 0.1	31	F 22 MAR 21 1.6	46	F 06 MAY 21 9.0	61		76	
2	F 22 MAR 21 wet	17	F 26 FEB 21 9.0	32	F 06 MAY 21 1.6	47	F 22 MAR 21 6.2	62		77	
3	B 05 APR 21 0.1	18	B 14 MAY 21 0.1	33	B 05 APR 21 9.0	48	C 19 APR 21 6.0	63		78	
4	C 19 APR 21 9.0	19	C 18 APR 21 TCT 4.0	34	B 14 MAY 21 6.0	49	B 03 MAY 21 3.0	64		79	
5	C 11 MAY 21 TCT 4.0	20	F 14 MAY 21 9.0	35	C 19 APR 21 0.1	50	F 16 APR 21 3.8	65		80	
6	F 08 FEB 21 0.1	21	F 22 MAR 21 5.0	36	F 14 MAY 21 6.2	51	F 06 MAY 21 8.0	66		81	
7	F 08 FEB 21 6.2	22	B 14 MAY 21 9.0	37	B 08 MAR 21 11.0	52	F 26 FEB 21 3.8	67		82	
8	C 11 MAY 21 MAC 3.0	23	F 06 MAY 21 0.1	38	B 08 MAR 21 0.1	53	F 24 MAY 21 wet	68		83	
9	B 16 SEP 20 11.0	24	F 22 MAR 21 wet	39	F 06 MAY 21 wet	54	C 11 MAY 21 CTS 5.0	69		84	
10	B 08 MAR 21 9.0	25	B 08 MAR 21 3.0	40	F 14 MAY 21 8.0	55	F 06 MAY 21 3.8	70		85	
11		26		41		56		71		86	
12		27		42		57		72		87	
13		28		43		58		73		88	
14	Duplicate	29	Duplicate	44	Duplicate	59	Duplicate	74	Duplicate	89	Duplicate
15	Spike	30	Spike	45	Spike	60	Spike	75	Spike	90	Spike

Remember to do a duplicate, a check standard, and a blank in each set of 10 samples.

Spiking Info:

Sample Info: